

HIDEX



Hidex **Sense** Gas mixer unit

425-387 Hidex Sense Gas mixer unit

Simultaneous regulation of CO₂ and O₂ environment enables hypoxia research, long-term cell proliferation measurements, and other biological studies of anaerobic organisms.

Hidex Sense environment control system consists of a high precision digital gas flow mixer connected to the enclosed Sense measurement chamber. The control unit is equipped with a touch screen to adjust settings and to monitor the gas concentrations.

The unit controls CO₂ and O₂ concentration by measuring and mixing CO₂ and Nitrogen continuously. Nitrogen is only needed for oxygen control. The actual gas flow and concentrations are automatically recorded by the mixer unit, from where these can be displayed on the instrument computer, or loaded to a USB device.



The system features special functionality to reduce gas consumption and fast recovery time after opening the microplate loading lid. An external audio alarm can also be connected to the unit in case the conditions require operator attention.

CO₂/O₂ Controller Setup

The controller regulates CO₂ and O₂ by infusing pure CO₂ and N₂ into the Sense measurement chamber, obtaining an adjustable CO₂ concentration in the range 0-20%, and an O₂ concentration in the range 1-20%.


The controller samples the CO₂ and O₂ gas concentrations every 30 second from the enclosed Sense microplate measurement chamber. The sampled gas is dehumidified with a membrane dryer (optional).



CO₂/O₂ Consumption and stability

The typical gas consumption is exceptionally low for a microplate reader, less than 0.5 liters/minute to keep CO₂ at a 5% steady state, and 7-10 liters/minute of Nitrogen to keep O₂ at a 1% steady state. Thus e.g. a 40 Kg CO₂ tank lifetime is estimated to 30 days of continuous use.

The graph shows typical CO₂ and O₂ concentration stability.



Sense CO₂/O₂ Controller Specifications

Operation mode	Adds pure CO ₂ and/or N ₂ to the measurement chamber
Concentration range	CO ₂ : 0-20%; O ₂ 1-20%
CO ₂ accuracy	0.1%
O ₂ accuracy	0.1%
Set point resolution	0.1%
Min/Max Pressure for CO ₂ input	700 mmHg to 5 barg (72.5 psig)
Min/Max pressure for N ₂ input	700 mmHg to 1.8 barg (26.1 psig)
Max infusion rate for CO ₂	12 l/min @ 5 barg (72.5 psig)
Max infusion rate for N ₂	27 l/min @ 1.8 barg (26.1 psig)
Outlet pressure	Ambient
CO ₂ Sensor	Non dispersive Infrared (NDIR) dual wave length detector
O ₂ Sensor	Optical sensor
Sensor life time	10 years for CO ₂ sensor, 5 years for O ₂
CO ₂ and O ₂ Input Gas connectors	6 mm OD push in fittings
Output Gas connector	Panel Mount Connector, 3.2 mm ID Tubing
Suction pump	Integrated pump for gas sampling
Gas sensor input connector	Panel Mount Connector, 3.2 mm ID Tubing
Dimensions	130x158x104 mm
Weight	1600 g
User interface	Touchscreen display
Sensor calibration mechanism	Comparison with external meter, or alternatively span gas
Recommended sensor calibration period	CO ₂ one year, O ₂ one year
Gas dryer (Optional)	To dehumidify the sampling gas prior to concentration measurement, needed only if the controlled volume is humid
Alarms	Buzzer, flashing and External Alarm Connector
Data logging	RS232 and MINI USB ports for data logging through Okolab DATA LOG software or any third party software



Hidex is a family owned high technology company which develops and manufactures high performance analysis equipment for life science research, radiation measurement and nuclear medicine. Our products utilize modern technology and excellent tradition of workmanship. With strong links to the scientific community we continue to innovate and develop to improve scientific research and safety of everyday life.

Today more than 3000 Hidex precision instruments are at service in leading laboratories worldwide as well as in some of the hardest conditions on the planet. Jungles and deserts, oil platforms and ocean going vessels – even submarines are no challenge for Hidex instruments.

Contact Hidex

Call us
Tel. +358 10 843 5570

Address
Lemminkäisenkatu 62
FIN-20520 Turku
Finland

E-mail
info@hidex.com
firstname.lastname@hidex.com

www.hidex.com